

# SETUP INSTRUCTIONS

## MyCareLink Smart™ Patient Monitor



### STEP 1

Have your smartphone or tablet readily available. Open the MyCareLink Smart Patient Monitor box and remove the contents:

- A Patient Reader
- 3 AA Batteries
- Setup Guide
- Patient Manual

The reader is designed to work with your Apple® or Android™ smartphone or tablet.



### STEP 2

Insert the 3 AA batteries into the back of the Reader by pushing the battery cover down and away from you. Make sure Bluetooth® within your smartphone or tablet is turned on.



### STEP 3

Using your smartphone or tablet, launch your browser and go to [mycarelinksmart.com](http://mycarelinksmart.com). Tap the download button, which will send you to either the Apple App store or the Google Play™ store—based on what type of smartphone or tablet you're using.



### STEP 4

Download the MyCareLink Smart App — just like any other app.



### STEP 5

Once the app is downloaded to your smartphone or tablet, open it and follow the prompts. Accept the terms and conditions; then, press the green button to begin.



### STEP 6

Turn the Reader on by pressing the gray button.

A continuous animation will be seen. That's because the app is trying to establish a Bluetooth connection with the Reader. You only need to press the gray button once.



## STEP 7

Locate the numeric code on the back of your Reader and enter it into the app using the app keypad. This action enables your Reader to work with your smartphone or tablet. You will also be prompted to use Bluetooth to pair your Reader with your smartphone or tablet. Simply accept this request.



## STEP 8

Place the Reader over your heart device. The Reader will display a green light and the app will display a green progress bar that lets you know the Reader is positioned correctly and receiving information from your device.




## STEP 9

The app will indicate your information is being sent to your clinic.



## STEP 10

A green check mark confirms that the transmission was successfully sent. Tap the home icon to go back to the home screen, which shows the date of your most recent transmission. At this point, you can exit the app. The Reader will turn off automatically.



**For more information, please consult the patient manual that was included with your MyCareLink Smart Monitor.**

Apple and the Apple logo are trademarks of Apple Inc., registered in the United States and other countries. App Store is a service mark of Apple Inc. Android and Google Play are trademarks of Google Inc. Bluetooth is a registered trademark of Bluetooth SIG, Inc.

Medtronic  
710 Medtronic Parkway  
Minneapolis, MN 55432-5604  
USA  
Tel: (763) 514-4000  
Fax: (763) 514-4879

Toll-free: 1 (800) 328-2518  
(24-hour technical support  
for physicians and medical  
professionals)

UC201602295 EN  
© Medtronic 2015.  
Minneapolis, MN.  
All Rights Reserved.  
Printed in USA. 09/2015

[medtronic.com](http://medtronic.com)

### Additional Information

The Medtronic MyCareLink Smart™ Reader is a prescription product indicated for use in the transfer of patient data from some Medtronic implantable cardiac devices based on physician instructions and as described in the product manual. The MyCareLink Smart Reader transmissions to the CareLink™ Network are subject to Wi-Fi and cellular service availability. This product is not a substitute for appropriate medical attention in the event of an emergency and should only be used as directed by a physician. The Medtronic CareLink Service is prescribed by your physician. This service is not for everyone. Please talk to your doctor to see if it is right for you. Your physician should discuss all potential benefits and risks with you. Although many patients benefit from the use of this service, results may vary. For further information, please call CareLink Patient Services at 1 (800) 929-4043 (8:00 a.m. to 5:00 p.m., Monday–Friday, Central time) or see the Medtronic website at [www.medtronic.com](http://www.medtronic.com).

**Medtronic**